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You may use the EL029 to connect a simple, single-pole switch (such as a micro-switch) to DrDAQ. You will need to obtain a switch and some insulated connecting wire. Remember to keep magnets or magnetic fields away from the EL029 when you use it with an external switch. Magnetic fields will still make the internal Reed Switch operate overriding the open setting of the external switch.

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Measure humidity

Setup the equipment

1. Clean the equipment. Make especially sure that the test container is clean; rinse it out with distilled water.
2. Put the test container into your insulated box.

Warning: The standard solutions may be harmful to your skin, eyes or when swallowed. Take all necessary precautions to avoid contact when preparing and using the standard solutions.

Caution: It is very important that the standard solution does not enter the case of the DD101. The solution may damage the electrical components of the sensor.

3. Put some of the standard solution into the test container. The solution should occupy about 5% of the volume of your container.
4. Put the stand into the test container. The stand should give you a clear platform above the level of the liquid. Do not allow any standard solution to spill on to the top of your stand.
5. Put the DD101 and DD100 sensors onto the stand. [If you are using a calibrated reference put this in too.]
6. Connect the DD101 and DD100 sensors to DrDAQ.

Make your measurements

7. Start the PicoLog Recorder software for DrDAQ. Make sure DrDAQ is receiving readings from the

sensors on Ext 1 and Ext 2.

8. Seal the test container and close the insulated box (if you are using a fan, start the fan).
9. Record the temperature and humidity inside the test container for at least one hour. You must wait for the temperature to stabilise and the DD101 to provide correct readings. This may take up to eight hours if you have used a large container.
10. Check that the DrDAQ plots for temperature and humidity have been constant over the last few minutes of your measurements.
[If you are using a calibrated reference make a reading of this now. Take care to minimise the changes in the setup so that the reading does not change significantly.]
11. Save your results.
12. Dispose of the solution as recommended by the supplier of the chemical.
Remember to take adequate precautions to protect your skin and eyes when disposing of the chemicals.
13. Do the measurements again for the other standard solutions. You should end up with a set of recorded measurements for each standard solution. Remember that you must provide at least two fixed points for calibration of the sensor.

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1. From your recorded results find the place near the end of the measurement time where the readings are most stable.
2. Write down these temperature and humidity readings for the standard solution in a table.
3. Look up the humidity that the standard solution should give for the temperature that you have recorded. Write this value in your table next to the value measured by the DD101. (The manufacturer of the standard solution should have provided you with a table for humidities for different temperatures).
4. Fill in entries for all the standard solutions you have measured.

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| Chemical | Measured humidity | Standard humidity | Temperature |
|--------------------|-------------------|-------------------|-------------|
| Potassium Chloride | 81% | 85% | 20 C |
| Magnesium Nitrate | 50% | 54% | 21 C |
| Magnesium Chloride | 30% | 33% | 20 C |

Keep this table safe for further reference. It will be useful if you have to reset the software with the PicoLog Recorder 'New Settings' command.

At this stage compare the measured and standard values. If you find that the measured values differ from the standard values by more than 10%, then there is something wrong. You may have had a non-saturated solution, some contamination in the solution, inadequate sealing of your test container, or possibly a damaged sensor. Check your calibration routine. If you still get large errors contact Pico for assistance.

[If you are using a calibrated reference this should give values very close to those quoted for your standard solutions. If these values disagree by more than a few percent suspect your calibration procedure. When the two are in agreement use the values given by the calibrated reference as the correct figures for your data table.]

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Open the PicoLog Parameter Scaling Dialog

1. Start PicoLog Recorder
2. Click on Settings
3. Click on Input Channels - the DrDAQ Measurements dialog appears
4. Select the input channel for the humidity sensor - for example 'External 1 Humidity'
5. Click on the Edit button - the Edit DrDAQ Measurements dialog appears

6. Click on the Options button - the Parameter Options dialog appears
 7. Click on the Scaling button - the Parameter Scaling dialog appears
- Enter your calibration data
8. Pull down the Scaling Method list and select 'Table lookup'
 9. To start your calibration table, click in the white rectangle - the text cursor appears
 10. Type in the value that you measured with the DD101 under 'Raw'. Type a space then enter the standard value under 'scaled'. Press return (the enter key).
 11. Enter all the pairs of measured and standard values that you have to build your table. You should end up with two columns of values.
Note: Your table must have at least two pairs of values.
- Click OK to close the dialogs
12. Click on the OK button to close the Parameter Scaling dialog - the Parameter Scaling dialog disappears
 13. Click on the OK button to close the Parameter Options dialog - the Parameter Options dialog disappears
 14. Click on the OK button to close the Edit DrDAQ Measurements dialog - the Edit DrDAQ Measurements dialog disappears
 15. Click on the OK button to close the DrDAQ Measurements dialog - the DrDAQ Measurements dialog disappears
- This completes the entry of the calibration data.

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Once you have successfully calibrated your DD101 write the date on the label and stick the label on the sensor (do not block the holes in the sensor with the label).

Testing Sensors

You can check whether a sensor is working properly or not by running through the calibration procedure. If you find that the sensor values differ from the standard values by more than 10% you may have a damaged sensor. In this case contact Pico for assistance.

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Removing Basic Probe Compensation

Open the PicoLog Parameter Scaling Dialog

1. Start the PicoLog Recorder software for DrDAQ. Make sure DrDAQ is receiving readings from the conductivity interface on Ext 1 (DD102).

2. Click on Settings
3. Click on Input Channels - the DrDAQ Measurements dialog appears
4. Select the input channel for the conductivity sensor - for example 'External 1 Conductivity'
5. Click on the Edit button - the Edit DrDAQ Measurements dialog appears
6. Click on the Options button - the Parameter Options dialog appears
7. Click on the Scaling button - the Parameter Scaling dialog appears

Removing the correction

8. Pull down the Scaling Method list and select 'None'
Note: You do not have to remove the calibration text from the box below.

Click OK to close the dialogs

9. Click on the OK button to close the Parameter Scaling dialog - the Parameter Scaling dialog disappears
10. Click on the OK button to close the Parameter Options dialog - the Parameter Options dialog disappears
11. Click on the OK button to close the Edit DrDAQ Measurements dialog - the Edit DrDAQ Measurements dialog disappears
12. Click on the OK button to close the DrDAQ Measurements dialog - the DrDAQ Measurements dialog disappears

This completes the removal of the probe correction factor.

2) Temperature Compensation

Conductivity in liquids varies with temperature. Expect a variation of around 2% of reading per °C difference from 25°C. Therefore at 23°C you would expect the conductivity to read around 4% lower than usually expected for that solution.

Overview

Temperature compensation is achieved using a DrDAQ DD100 temperature sensor connected to your spare EXT input. PicoLog is configured to take readings from both the DD102 and DD100 sensors. A calculated parameter is then created in PicoLog which allows the temperature readings from DD100 to adjust the readings produced by DD102.

Procedure

1. Start the PicoLog Recorder software for DrDAQ. Make sure DrDAQ is receiving readings from the sensors on Ext 1 (DD102) and Ext 2 (DD100).
2. Perform the "Basic Probe Compensation" described above.

Create a calculated parameter

3. Click on Settings
4. Click on Calculated Parameters - the Calculated Parameters dialog appears
5. Click on the Add button - the Edit Calculated Parameters dialog appears
6. Type in a name for the calculated parameter - for example 'Temp Compensated Conductivity'
7. Pull down the Input Parameter box 'A' and then select 'External 1 Conductivity'
8. Pull down the Input Parameter box 'B' and then select 'External 2 Temperature'
9. Click in the Equation box - the text cursor appears
10. Type in the following equation $A+(A*0.02*(25-B))$ (see below for explanation)
11. Click on the Options button - the Parameter Option dialog appears
12. In the Units box type uS
13. Click in the Decimal Places box and set to zero.

Click OK to close the dialogs

14. Click on the OK button to close the Parameter Scaling dialog - the Parameter Scaling dialog disappears
15. Click on the OK button to close the Edit Calculated Parameters dialog - the Edit Calculated Parameters dialog disappears
16. Click on the OK button to close the Calculated Parameters dialog - the Calculated Parameters dialog disappears

This completes the temperature compensation.

Explanation of above equation :

- A is the raw conductivity reading from DD102

- B is the temperature reading from DD100
- 0.02 equates to 2% per °C variation in conductivity with temp - this may be different for your liquid.
- 25 is the reference temperature in °C - this may be different for your liquid.

3) Full Calibration

Basis of the calibration method

Calibration relies on the general physical properties of calibration solutions. The conductivity in such solutions is known quite accurately and is used to calibrate the sensors. Solutions for this purpose are readily available.

Note: that the conductivity that exists in all such solutions varies with temperature so accurate temperature compensation must be used.

Warning: The standard calibration solutions may be harmful to your skin, eyes or when swallowed. Take all necessary precautions to avoid contact when preparing and using the standard solutions.

Note: For the highest accuracy you should compare the readings from your DD102 with the readings from a calibrated 'laboratory standard' reference conductivity gauge.

Calibration of the sensor involves:

- Preparation of the calibration equipment.
- Measurement of the standard calibration solutions and temperature.
- Creation of calibration data for the DrDAQ software.
- A check that the calibrated sensor is accurate.

Note: You must measure at least two different standard solutions to provide two or more fixed points for calibration.

Equipment you will need

- DrDAQ with DD102 sensor
- Pico DD100 Temperature Sensor
- At least two Standard Solutions
- An insulated box
- A label for the DD102

[Optional: high accuracy, calibrated reference 'laboratory standard' conductivity gauge]

Warning: It is very important that the standard solution does not enter the case of the DD102. The solution may damage the internal electrical components.

Standard Solutions

Whichever solutions you choose make sure that the conductivity they give is within the operating range of the DD102 sensor (that is: greater than 70uS and less than 20000uS conductivity). You should use at least two standard solutions to give two fixed points for your calibration. Standard reference books such as Kaye & Laby "Tables of Physical and Chemical Constants" (Longman) give tables with the conductivity of standard solutions.

Insulated box

Because the conductivity produced by the standard solution depends on temperature, you should consider putting your test solutions inside an insulated box to ensure the temperature of the solution is stable during the calibration procedure. A domestic cool box (without the cold blocks) is suitable.

Preparation

Prepare your standard solutions before you start the calibration procedure. Make sure that all the containers you use for preparation are thoroughly clean before use, contamination of the solution will alter the conductivity. Allow time for the solution to stabilise at room temperature before use. Always use fresh

solutions to ensure that the chemicals have not become contaminated or degraded.

Allow the standard solutions and other equipment to reach a stable temperature (ideally 25 degrees Celsius) before you start.

Measure conductivity

Setup the equipment

1. Clean the equipment. Make especially sure that the solution containers are clean; rinse out with distilled water if required.
2. Put the solution containers into your insulated box.
3. Connect the DD102 and DD100 sensors to DrDAQ.
4. Start the PicoLog Recorder software for DrDAQ. Make sure DrDAQ is receiving readings from the sensors on Ext 1 and Ext 2.
5. Perform the "Basic Probe Compensation" described above.
6. Perform the "Temperature Compensation" described above.
7. Put the conductivity probe into the first solution. [If you are using a calibrated reference put this in too.]

Make your measurements

8. Monitor the conductivity and temperature of the solution for at least one minute while agitating the probes. You must wait for the readings to stabilise before taking your calibration measurement.
9. Once the conductivity and temperature readings have stabilised, take a reading for both
10. [If you are using a calibrated conductivity reference make a reading of this now. Take care to minimise the changes in the setup so that the reading does not change significantly.]
11. Save your results.
12. Dispose of the solution as recommended in the solution instructions. Remember to take adequate precautions to protect your skin and eyes when disposing of the chemicals.
13. Repeat the measurements for the other standard solutions. You should end up with a set of recorded measurements for each standard solution. Remember that you must provide at least two fixed points for calibration of the sensor.

Create calibration data

When the measurements are completed you need to make a written calibration table.

1. Enter the temperature compensated conductivity that you have measured in the "Measured Conductivity" column.
2. Enter the conductivity of the solution you have tested in the "Solution Conductivity" column.
3. Fill in entries for all the standard solutions you have measured.

You should end up with a table something like this:

| Type of Chemical Solution | Measured Conductivity | Solution Conductivity |
|----------------------------|-----------------------|-----------------------|
| Low Conductivity Solution | 89uS | 84uS |
| Med Conductivity Solution | 1350uS | 1413uS |
| High Conductivity Solution | 11573uS | 12880uS |

Keep this table safe for future reference. It will be useful if you have to reset the software with the PicoLog Recorder 'New Settings' command.

At this stage compare the measured and expected values. If you find that the measured values differ from the standard values by more than 10%, then there is something wrong. You may have a contaminated solution, an incorrect temperature reading, a damaged probe or possibly a damaged interface. Check your calibration routine. If you still get large errors contact Pico for assistance.

[If you are using a calibrated reference sensor, this should give values very close to those quoted for your standard solutions. If these values disagree by more than a few percent suspect your calibration procedure. When the two are in agreement use the values given by the calibrated reference as the correct figures for your data table.]

Once you have your calibration data you must enter it into the Pico software to calibrate the sensor:

Edit the calculated parameter

1. Click on Settings
2. Click on Calculated Parameters - the Calculated Parameters dialog appears
3. Click on the calculated parameter you entered earlier - called 'Temp Compensated Conductivity'
4. Click on the Edit button - the Edit Calculated Parameters dialog appears
5. Click on the Options button - the Parameter Option dialog appears
6. Click on the Scaling button - the Parameter Scaling dialog appears

Enter your calibration data

7. Pull down the Scaling Method list and select 'Table lookup'
8. To start your calibration table, click in the white rectangle - the text cursor appears
9. Using the table you created above, type in the "measured conductivity" under 'Raw'. Type a space then type the "solution conductivity" value under 'scaled'. Press return (the enter key).
10. Enter all the pairs of measured and solution values that you have to build your table. You should end up with two columns of values.

Note: Your table must have at least two pairs of values.

Click OK to close the dialogs

11. Click on the OK button to close the Parameter Scaling dialog - the Parameter Scaling dialog disappears
12. Click on the OK button to close the Parameter Options dialog - the Parameter Options dialog disappears
13. Click on the OK button to close the Edit Calculated Parameters dialog - the Edit Calculated Parameters dialog disappears
14. Click on the OK button to close the Calculated Parameters dialog - the Calculated Parameters dialog disappears

This completes the entry of the calibration data.

Note: You can remove the calibration data by opening the Parameter Scaling dialog and setting the pull down list to 'none'. If you want to completely delete the calibration scaling you can either edit and delete the entries in the Parameter Scaling dialog or use the 'New Settings' command from the main menu. If you use the '**New Settings**' command **all scaling data will be lost**.

Check the calibration

To be sure that the calibration has been successful you must repeat the measurement stage of the procedure. When the check measurements have been completed there should be very close agreement between the measured and standard values (that is within 5% of the value). If this is not the case check that you have entered the calibration data correctly and repeat the process.

Once you have successfully calibrated your DD102 write the date on the label and stick the label on the sensor. The output of the sensor may vary over time. For this reason you should calibrate the sensor regularly.

Testing the interface and probe

You can check whether the interface and probe combination is working properly by running through the calibration procedure. If you find that the sensor values differ from the standard values by more than 10% you may have a damaged probe or even a damaged interface. In this case contact Pico for assistance.

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Each external sensor socket has two channels. One is used to automatically detect the type of sensor: the other is an analogue input that represents the sensor reading.

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The DrDAQ driver has built-in scaling for each of the built-in and pico-supplied sensors. You can incorporate scaling for your own sensors by adding a file called `drdaq.dds`, which contains details of your sensor.

The values returned by the driver are integers that represent fixed-point decimal number. For example, the driver treats pH as a value with two decimal places, so a pH of 7.65 is returned as 765.

You can call the routine `drdaq_get_channel_info` to find out how many decimal places a channel is using, and also to get a divider that converts the integer value to the corresponding real number. For pH, the returned divider is 100, so 765 divided by 100 gives 7.65.

For some sensors, there is more than one possible scaling available. You can call `drdaq_get_scalings` to get a list of valid scaling codes, then call `drdaq_set_scaling` to select one of them. Once selected, `drdaq_get_channel_text` and `drdaq_get_channel_info` will return full information about the selected scaling. If you do not use `drdaq_set_scaling`, the driver will automatically select the first available scaling for each channel.

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Channels

The following table gives details of the scaling for each channel.

| Channel number | Input | Min Value | Max Value | Decimal Places | Units |
|----------------|----------------|-------------------|-----------|----------------|----------------|
| 1 | Sound waveform | -100 | 100 | 1 | |
| 2 | Sound level | 50 | 100 | 2 | dBa |
| 3 | Voltage | 0 | 5000 | 0 | mV |
| 4 | Resistance | 0 | 1000 | 1 | kOhm |
| 5 | pH | 0 | 14 | 2 | |
| 6 | Temperature | 0 | 100 | 1 | C ² |
| 7 | Light | 0 | 100 | 1 | |
| 8 | External 1 | Depends on sensor | | | |
| 9 | External 2 | | | | |

¹ Fahrenheit and Kelvin are also available- see drdaq_apply_fix

² Fahrenheit and Kelvin are also available- see drdaq_apply_fix

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drdaq_set_unit

```
PREF1 short PREF2 drdaq_set_unit (short port);
```

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drdaq_set_led

```
PREF 1 void PREF2 drdaq_set_led (short do_value);
```

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drdaq_get_value

PREF 1 short PREF2 drdaq_get_value (short channel);

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drdaq_get_value_and_time

```
PREF 1 void PREF2 drdaq_get_value_and_time (  
    short channel,  
    unsigned long * sample_time,  
    short * value);
```

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drdaq_set_trigger

```
PREF1 void PREF2 drdaq_set_trigger ( short enabled,  
                                   short auto_trigger,  
                                   short auto_ms,  
                                   short channel,  
                                   short dir,  
                                   short threshold,  
                                   short delay);
```

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drdaq_set_interval

```
PREF1 unsigned long PREF2 drdaq_set_interval (  
    unsigned long us_for_block,  
    unsigned long ideal_no_of_samples,  
    short * channels,  
    short no_of_channels);
```

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```
int channels [3];  
channels [0] = 2;  
channels [1] = 3;  
channels [2] = 5;
```

```
drdaq_set_interval (10000, 100, channels, 3);
```

The routine returns the actual time to collect this number of samples. This actual time may be greater than the target time if you specified a sampling interval that is faster than your computer can manage. If the specified sampling rate was too fast, you have the following choices:

- if the total time is important, collect fewer than the ideal number of samples so that the total block time is correct

- if the number of samples is important, collect the same number of samples then allow for the fact that they took longer to collect.

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drdaq_get_times_and_values

```
PREF1 unsigned long PREF2 drdaq_get_times_and_values (  
    long HUGE * times,  
    short HUGE * values,  
    unsigned long no_of_values);
```

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drdaq_get_unit_info

```
PREF1 short PREF2 drdaq_get_unit_info (  
    char * str,  
    short str_lth,  
    short line,  
    short port);
```

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drdaq_get_channel_text

```
PREF1 short PREF2 drdaq_get_channel_text (  
    char * str,  
    short channel,  
    short field);
```

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drdaq_get_channel_info

```
PREF1 short PREF2 drdaq_get_channel_info (  
    short * min_value,  
    short * max_value,  
    short * places,  
    short * divider,  
    short * is_fast,  
    short channel);
```

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drdaq_get_scalings

```
PREF1 short PREF2 drdaq_get_scalings (  
    short * ids,  
    short channel);
```

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drdaq_set_scaling

```
PREF1 short PREF2 drdaq_set_scaling (  
    short channel,  
    short scaling);
```

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